



Modul-Bio

High-throughput solutions

- BARCODE SOLUTIONS**
- BARCODE LABELS**
- BARCODE SOFTWARE**
- BARCODE PRINTERS**
- BARCODE SCANNERS**

Identification solutions for long term storage of samples in liquid nitrogen

- Technology used on four continents by repositories, government agencies, universities or private sector
- Millions of tubes labelled for storage in countless liquid nitrogen tanks
- Labels proven effective for over 6 years


BARCODE LABELS for long term storage of samples in liquid nitrogen



Used by the Center for Disease Control (C.D.C.)
Atlanta, GA - USA.


- Secure "wrap-around" design protects barcode
- Label already frozen vials
- Effective at extreme temperatures (-196° to +90°C)
- Solvant resistant
- Available in multiple sizes

white area
for printing



transparent area
for content
visualization

BARCODE SOFTWARE



- Design and print your own labels
- Information on labels may include :
 - Barcode and human readable
 - Type of samples (DNA, RNA, Blood, Tissue, ...)
 - Date, Aliquot n#, Laboratory name, etc.
- Available with database access

BARCODE PRINTERS - ZEBRA™

- Model Z4M+, 105SL, 110XiIII+
- High resolution 300 dpi and higher
- Network card



BARCODE SCANNERS - SYMBOL™

- Model LS 2208, LS 4004
- High resolution scanners



Contact : SAS Modul-Bio
232, bd Sainte-Marguerite - 13009 Marseille - FRANCE
Voice : + 33 (0)4 91 22 33 05 - Fax : + 33 (0)4 91 22 35 72
info@modul-bio.com - www.modul-bio.com

Contact us for more information

