

PUBLISH & SHARE YOUR SAMPLE COLLECTIONS ON THE WEB



PERSONALISED ONLINE BIOSPECIMEN CATALOGUE

A software solution for centralised sharing of biospecimens from multiple biobanks

- Presentation of biobanks sharing biospecimen.
- Each biobank can securely upload a file containing biospecimen information.
- Powerful search function for selecting biospecimens.
- Management of sample requests & shipping.
- Integration of sample results from researchers.

Visibility

Better visibility for your biobank.

Distribute more biospecimens


Increase collaborations with researchers.

Interoperability

Publish sample data from different datasources.

Personalised

Can be adapted to your own design (CSS).



More than 40 biobanks use eMBioBANK®

BIOBANKS - HOSPITALS - RESEARCH CENTERS - BIOTECHNOLOGY & PHARMACEUTICAL COMPANIES

↘ eMBioBANK® features



Biobank presentation

Presentations published for each biobank



Share biospecimens

Present multiple collections in a single online catalogue



User profiles

Create a user profile for each biobank
Allows flexible role & privilege access



Sample data transfer

Secure uploading of data
Quality control checks
Sample & Patient anonymisation
Compatible with MBioLIMS®



Searches

Complex & multi-criteria searches
Save custom searches



Sample Request Orders

Researchers can submit sample requests
Validation of requests
Researchers can track request processing



Sample Shipping

Record samples shipments
Generates packing slips



Security

HTTPS protocol, sample anonymisation and secure data file transfer



Record data from researchers

Save a variety of data (genomic, proteomic, quality control, images, files, etc.) returned by researchers



Reporting

Powerful reporting and statistical tools



International standards

Classifications included: ICD10, CIM-O, SNOMED, etc.

↘ Technical specifications:

- Follows OECD guidelines, FDA 21 CFR part 11, GLP.
- Robust (Web2.0, Java, Oracle).
- 100% web-based software with secured access (HTTPS).
- Works with all browsers (Explorer, Firefox, Chrome, Safari) and with Linux, Windows, OS X.
- Separate server solution.
- Hosting solutions.

Discover our solutions

MBioLABEL

MBioLIMS

eMBioBANK



Modul-Bio
OPTIMIZE YOUR SAMPLES

Contact us for a demonstration:
demo@modul-bio.com

www.modul-bio.com